
Delto Gripper

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Ver. 1.0.0

INFORMATION

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4.1 URCap Installation

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DG-5F

To use the URCaps feature of the Delto Gripper-5F, you first need to connect the gripper.

1. Set the left switch to ① and the right switch to ③ RS485 or ④ TCP/IP before powering the gripper.

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
DG-5F

To use the URCaps feature of the Delto Gripper-5F, you first need to connect the gripper.

3. Change the value and test it.


1 **DISCONNECT**

2 **Basic** **Advanced** **etc**

3 Initial Pose Recipe Number:
 Initial Gain Recipe Number:
 Initial Grasp Recipe Number:

4 **APPLY** **WRITE** 5 PLY: Applies the values to the gripper. If WRITE is not used, the applied values will be reset upon reboot.
 WRITE: Saves the values applied to the gripper to the EEPROM. The values will be retained even after a reboot.

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4.1.2.1 Basic

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	Basic	Advanced	etc
1	Initial Pose Recipe Number:	1	
2	Initial Gain Recipe Number:	1	
3	Initial Grasp Recipe Number:	1	
	APPLY WRITE		4 SYSTEM OFF

* **APPLY**: Applies the values to the gripper. If **WRITE** is not used, the applied values will be reset upon reboot.
 * **WRITE**: Saves the values applied to the gripper to the EEPROM. The values will be retained even after a reboot.

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4.1.2.2 Advanced

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	Basic	Advanced	etc
1	Settings Joint Offset Joint Offset Joint Inpose Moving Inpose Slave ID IP Address Subnet Mask Gateway Port Low Pass Filter	Joint 1 0.0	0.0

* **APPLY**: Applies the values to the gripper. If **WRITE** is not used, the applied values will be reset upon reboot.
 * **WRITE**: Saves the values applied to the gripper to the EEPROM. The values will be retained even after a reboot.

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4.1.2.3 etc.

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	Basic	Advanced	etc
1	Product ID : 5F12		
2	Firmware Version : 2.2		

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4.2.1 CURRENT

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DG-5F CURRENT

Set the gripper to current-based position control mode.

1
☒ Activate
 ☐ Deactivate

2

Joint 1 (mA)	Joint 2 (mA)	Joint 3 (mA)	Joint 4 (mA)
0	0	0	0
Joint 5 (mA)	Joint 6 (mA)	Joint 7 (mA)	Joint 8 (mA)
0	0	0	0
Joint 9 (mA)	Joint 10 (mA)	Joint 11 (mA)	Joint 12 (mA)
0	0	0	0
Joint 13 (mA)	Joint 14 (mA)	Joint 15 (mA)	Joint 16 (mA)
0	0	0	0
Joint 17 (mA)	Joint 18 (mA)	Joint 19 (mA)	Joint 20 (mA)
0	0	0	0

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
4.2.2 GAIN

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DG-5F GAIN

The screenshot shows the 'DG-5F GAIN' configuration interface. It includes a control mode selector (1) with radio buttons for 'PD Control' (selected) and 'PID Control'. A dropdown menu (2) is set to 'Finger 1' and 'P Gain'. A large input area (3) for joint gains (P, D, I, L) for joints 1 through 4, with '1' entered for all P gains and '0' for all D, I, and L values. A checkbox (4) for 'Advanced Settings' is checked. A 'Load Gain' section (5) has 'ACTIVE' and 'DEFAULT' buttons. A 'Save Recipe' section (6) has an empty input field and a 'SAVE' button. An 'EEPROM Write' section (7) has a 'WRITE' button. A 'TEST' button (8) is at the bottom right. Red warning text is present below the 'SAVE' and 'WRITE' buttons.

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4.2.3 GPIO

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DG-5F GPIO

This node configures the GPIO operation of the gripper.
It allows you to set the value of a GPIO output or check the status of a GPIO input.

1

☒ Set GPIO Output
☐ Output 1 ☐ Output 2 ☐ Output 3

2

TEST

3

☐ Wait for GPIO Input LOW ▼

4

☐ Wait time before node start is 0.0 seconds
☐ Wait time after node completion is 0.0 seconds

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3

4 Position Control

Finger 1	Finger 2	Finger 3	Finger 4	Finger 5
☐ Joint 1	☐ Joint 5	☐ Joint 9	☐ Joint 13	☐ Joint 17
☐ Joint 2	☐ Joint 6	☐ Joint 10	☐ Joint 14	☐ Joint 18
☐ Joint 3	☐ Joint 7	☐ Joint 11	☐ Joint 15	☐ Joint 19
☐ Joint 4	☐ Joint 8	☐ Joint 12	☐ Joint 16	☐ Joint 20

5 ☒ Update Grasp Option

6 ☒ Advanced Settings

7 Save Recipe

SAVE

8 EEPROM Write

WRITE

* Warning! The existing saved data will be overwritten. Please use with caution.

* To use the saved recipe after a reboot, please perform **[WRITE]**.

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
4.2.5 JOINT

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DG-5F JOINT

The screenshot shows the DG-5F JOINT control interface. It features a grid of 20 joint value inputs (Joint 1 to Joint 20), each with a numeric field and a range (0.0 / 300). Below this is a section for unit selection (Degree(*) / Time(ms)) and a checkbox for Advanced Settings. The interface is divided into two main sections: 'Load Joint' and 'Teaching Mode'. The 'Load Joint' section includes buttons for 'TARGET' and 'CURRENT'. The 'Teaching Mode' section includes buttons for 'ACTIVATE' and 'DEACTIVATE'. Below these are sections for 'Save Recipe' (with 'TARGET' and 'CURRENT' buttons) and 'EEPROM Write' (with a 'WRITE' button). A 'TEST' button is located at the bottom right. Numbered callouts 1 through 8 highlight specific elements: 1 points to the joint value inputs, 2 points to the unit selection, 3 points to the Advanced Settings checkbox, 4 points to the 'TARGET' button, 5 points to the 'DEACTIVATE' button, 6 points to the 'CURRENT' button, 7 points to the 'WRITE' button, and 8 points to the 'TEST' button.

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
4.2.6 LOAD

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DG-5F LOAD

1 ☒ Recipe Pose ☐ Recipe Gain ☐ Recipe Grasp ☐ Recipe Blend

2 Load Recipe Number: 1

3 ☒ Advanced Settings

4 Save Blend Recipe
Press [CLEAR] to initialize

Blend Pose Recipe Number: 1

Blend Gain Recipe Number: 1

Blend Wait Time: 100

5 CLEAR 6 ADD

7 Blend Index: SAVE

* Warning! The existing saved data will be overwritten.
Please use with caution.

EEPROM Write

8 WRITE

To use the saved blend recipe after a reboot,
please perform [WRITE].

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4.2.7 WAIT

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DG-5F WAIT

Select the gripper state condition that will trigger the program's next action.

1

☒ Wait until the joint stops moving
☐ Wait until the joint reaches the target value
☐ Wait until the gripper's Blend function becomes STOP ▼

2

☐ Wait time before node start is 0.0 seconds
☐ Wait time after node completion is 0.0 seconds

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© 2025 주식회사 테솔로 All rights reserved